




SPECIAL PROMOTIONAL OFFERS: May 1 – August 31, 2019

Post-it® Custom Printed BIG Pads		Item #	Imprint	LOW 25 QTY!		50-99	100-249	250-499			
8" x 8" 20 self-stick sheets		PG88	1 spot color**	9.15	5.92	4.99	3.86	3.19	2.84	2.39	c
11-3/4" x 11-3/4" 20 self-stick sheets		PG1175	1 spot color**	15.25	11.07	9.85	8.19	6.99	6.15	5.10	c
15-3/4" x 15-3/4" 20 self-stick sheets		PG1575	1 spot color**	16.75	14.37	12.25	9.73	7.79	7.16	5.79	c

MADE IN THE USA
WITH GLOBALLY SOURCED MATERIALS

Be sure to reference offer **#P4277** on your P.O.

SAVE ON **BIG IDEAS** WITH **BIG PADS!**



7 Business Days Production

Available in limited paper colors. Call for details. PMS Ink Color Match at no additional charge.** Black ink recommended for Neon/Ultra papers.

U.S. Products and Pricing. Quantity price listed per unit. Standard catalog upcharges apply.
No additional charge for Neon/Ultra and Recycled papers. Two (2) adhesive stripes mandatory on all sizes. Super sticky and two-sided printing not available.
**No charge for PMS inks except fluorescent, special mix and metallic



Spring/
Summer



2019

Special Offers Valid May 1 – August 31, 2019
3M Promotional Products


asi 91240 | ppai 113638 | SAGE® 55850 | UPIC POST-IT | pppc 840606
3M, Post-it and the color Canary Yellow are trademarks of 3M. © 3M 2019. Printed in U.S.A.



Special Promotional Offers

May 1 – August 31, 2019

1-C SHEET
IMPRINT
INCLUDED!



Post-it® Custom Printed Notes Cubes		Imprint	250	500	1000	
C345 Half Cube 2-3/4" x 2-3/4" x 1-3/8"		1 spot color**, 1 design side print	3.38	2.89	2.71	c
		2 spot colors**, 1 design side print	3.49	2.99	2.78	
		4-color process, different design each side (4 designs total!)*	5.19	4.49	3.90	
C690 Full Cube 2-3/4" x 2-3/4" x 2-3/4"		1 spot color**, 1 design side print	5.31	4.92	4.59	c
		2 spot colors**, 1 design side print	5.45	4.97	4.66	
		4-color process, different design each side (4 designs total!)*	6.89	6.29	5.69	
1 spot color**, individual sheet print (no bleeds)			INCLUDED ALL QUANTITIES			
9 Business Days Production		Standard White paper only. Other paper options available, call for details.				

MADE IN THE USA
WITH GLOBALLY SOURCED MATERIALS

Be sure to reference offer **#P4277** on your P.O.

MATCH YOUR EXACT LOGO PMS COLOR!



Post-it® Custom Printed Notes		Item #	Sheet Count	Imprint	500-999	1000-1999	2000-2499	2500-4999	
3 x 4		32	25-sheets	1 spot color**	1.79 .58	.59 .42	.54 .39	.40 .33	c
				2 spot colors**	1.26 .80	.79 .52	.73 .49	.43 .37	
	35	50-sheets	1 spot color**	1.27 .84	1.00 .68	.93 .63	.65 .52	c	
			2 spot colors**	1.73 1.01	1.23 .81	1.13 .74	.69 .57		
4 x 6		42	25-sheets	1 spot color**	1.23 .87	.99 .67	.91 .64	.64 .52	c
				2 spot colors**	1.76 1.12	1.19 .84	1.10 .76	.70 .59	
	45	50-sheets	1 spot color**	1.92 1.26	1.61 1.01	1.48 .95	1.07 .85	c	
			2 spot colors**	2.29 1.49	1.73 1.15	1.61 1.07	1.12 .91		
7 Business Days Production				PMS Ink Color Match at No Additional Charge**					

MADE IN THE USA
WITH GLOBALLY SOURCED MATERIALS

Be sure to reference offer **#P4277** on your P.O.

OUR BEST VALUE FOR NOTES!



Post-it® Custom Printed Notes		Item #	Sheet Count	Imprint	LOW 250 QTY!	500	1000	2500	
2-3/4 x 3		PD33P-25	25-sheets	4-color process*	.80	.62 .52	.39 .35	.35 .31	c
		PD33P-50	50-sheets		.91	.89 .65	.66 .49	.61 .45	
3 x 4		PD34P-25	25-sheets	4-color process*	.82	.71 .56	.41 .35	.37 .32	c
		PD34P-50	50-sheets		1.02	.97 .79	.73 .58	.67 .54	
3 x 5		PD35P-25	25-sheets	4-color process*	1.17	.80 .68	.57 .48	.52 .44	c
		PD35P-50	50-sheets		1.59	1.31 1.04	.90 .79	.87 .69	
4 x 6		PD46P-25	25-sheets	4-color process*	1.33	1.05 .88	.73 .66	.67 .61	c
		PD46P-50	50-sheets		1.74	1.65 1.24	1.26 .97	1.08 .94	
5 Business Days Production				Standard White paper only. Orders must be placed in increments of 250.					

MADE IN THE USA
WITH GLOBALLY SOURCED MATERIALS

Be sure to reference offer **#P4278** on your P.O.

24 GET THEM SOONER!
ALL sizes of full-color notes!
By request. Details below.

Low Quantity Design Change Program Receive TOTAL QUANTITY pricing on design changes AS LOW AS 50 NOTE PADS PER DESIGN. Add **\$16.50 (G)** per design change. Minimum order quantity is 250 total note pads. Design changes must be ordered in increments of 50 note pads. Full-Color printing only.

U.S. Products and Pricing. Quantity price listed per unit. Standard catalog upcharges apply.

*For Full-color printing, requests calling for a specified PMS color will be converted to the closest possible 4-color process match. Exact color match cannot be guaranteed.
**No additional charge for PMS inks except fluorescent, special mix and non-standard metallic (\$285.00 (G) per color). Not available for all products.

24 24-Hour Service (available by request):
1. 24-Hour Service must be in subject line of email to 24hour@3Mpromote.com. **2.** 24-Hour Service must be clearly stated on Purchase Order. **3.** Orders may be faxed to 1-800-625-2685. **4.** 3M must receive a complete Purchase Order and production-ready artwork by 1:00PM CST, your order ships the next business day. **5.** Missing or conflicting information or an order requiring artwork modifications may delay the order. **6.** If a proof is required for a 24-Hour Service order, service starts after the proof is approved. **7.** Weekend and holidays are not considered business days. **8.** 24-Hour Service is defined as the time it takes 3M to print your order after you have approved your proof and does not include the time an order takes to be delivered once it leaves our production facility. Be sure to account for delivery times when you place your order. **9.** 3M cannot guarantee 24-Hour Service in the event of severe weather, natural disasters, or any other events outside of 3M's reasonable control. **10.** Must be an active account with good credit or have valid credit card on file. **11.** Products available for 24-Hour Service: PD33P-25, PD33P-50, PD34P-25, PD34P-50, PD35P-25, PD35P-50, PD46P-25 and PD46P-50. **12.** Drop ships, design change and copy change not available.

NO ADDITIONAL CHARGE for Virtual Samples • Set-up Ink Screens • PMS Inks**